




























PMI Non-Stick PTFE Coated Electrodes

PMI disposable electrodes resist splitting in high temperatures and minimize the probability of damage to surrounding tissue. PTFE insulation withstands flaking at high temperatures and resists eschar buildup during surgical procedures. PMI cost-saving electrodes are available in several configurations to assist in various electrosurgical procedures.

ELECTROSURGERY TIPS

-  P0009 5MM BALL, 5" (12.70cm)
-  P0012 2.50 (6.40cm) BLADE
-  P0012A 2.75" (6.99cm) BLADE
-  P0012AM 2.75" (6.99cm) MODIFIED BLADE
-  P0012M 2.50" (6.99cm) MODIFIED BLADE
-  P0013 2.75" (6.99cm) NEEDLE
-  P0013M 2.75" (6.99cm) MODIFIED NEEDLE
-  P0014 6.5" (16.51cm) BLADE
-  P0014A 4" (10.16cm) BLADE
-  P0014AM 4" (10.16cm) MODIFIED BLADE
-  P0014M 6.5" (16.51cm) MODIFIED BLADE
-  P0015 5MM BALL, 2" (5.08cm)
-  P0016A 4" (10.16cm) NEEDLE
-  P0016AM 4" (10.16cm) MODIFIED NEEDLE
-  P0016M 6.5" (16.51cm) MODIFIED NEEDLE

LAPAROSCOPIC

-  P0017 LAPAROSCOPIC STD BLADE 33cm
-  P0018 LAPAROSCOPIC SPATULA 33cm
-  P0019 LAPAROSCOPIC CURVED BLADE 33cm
-  P0020 LAPAROSCOPIC L-HOOK 33cm
-  P0021 LAPAROSCOPIC J-HOOK 33cm
-  P0022 LAPAROSCOPIC NEEDLE 33cm
-  P0024 LAPAROSCOPIC BALL 33cm
-  P0027 LAPAROSCOPIC THIN CURV SPATULA 33cm
-  P0100 LAPAROSCOPIC WIRE L-HOOK 33cm
-  P0105S LAPAROSCOPIC WIRE J-HOOK SPLIT 33cm


Split-Stem Configurations Available








Item numbers ending in **S** are offered in a *split-stem* configuration for use with reusable foot-control cords.

PUSH BUTTON PENCIL W/ COATED BLADE & HOLSTER

P0035H



LLETZ

-  P0410
10mm X 7mm
LLETZ SQUARE 13cm
-  P0450
10mm X 10mm
LLETZ LOOP 13cm
-  P0460
15mm X 12mm
LLETZ LOOP 13cm
-  P0470
20mm X 12mm
LLETZ LOOP 13cm
-  P0480
20mm X 15mm
LLETZ LOOP 13cm

TUNGSTEN MICRO NEEDLE

-  P1650 MICRO NEEDLE 3cm S1 RAI GHI
-  P1651 MICRO NEEDLE 4cm STRAIGHT
-  P1652 MICRO NEEDLE 3cm 45° ANG,
3mm FROM TIP
-  P1653 MICRO NEEDLE 3cm 45° ANG,
1cm FROM TIP

SINGLE USE REM GROUNDING PAD W/CORD

P4500



Item #	Description	Qty.
STANDARD COATED, GREY		
P0009	5MM BALL, 5" (12.70cm)	12
P0012	2.50(6.40cm) BLADE	12
P0012A	2.75" (6.99cm) BLADE	12
P0012AM	2.75" (6.99cm) MODIFIED BLADE	12
P0012M	2.50 (6.99cm) MODIFIED BLADE	12
P0013	2.75" (6.99cm) NEEDLE	12
P0013M	2.75" (6.99cm) MODIFIED BLADE	12
P0014	6.5" (16.51cm) BLADE	12
P0014A	4" (10.16cm) BLADE	12
P0014AM	4" (10.16cm) MODIFIED BLADE	12
P0014M	6.5" (16.51cm) MODIFIED BLADE	12
P0015	5MM BALL, 2" (5.08cm)	12
P0016AM	4" (10.16cm) MODIFIED NEEDLE	12
P0016A	4" (10.16cm) NEEDLE	12
P0016M	6.5" (16.51cm) MODIFIED NEEDLE	12
P0028	5.80" BAYONET NEEDLE 14.73cm	12
P0028M	5.80" MODIFIED BAYONET NDL 14.73cm	12
P0113	2.75" (6.99cm) BLUNT NEEDLE	12
P0113M	2.75" (6.99cm) BLUNT NEEDLE MOD	12
TUNGSTEN MICRO NEEDLE UNCOATED		
P1650	MICRO NEEDLE 3cm STRAIGHT	10
P1651	MICRO NEEDLE 4cm STRAIGHT	10
P1652	MICRO NEEDLE 3cm 45° ANG, 3mm from TIP	10
P1653	MICRO NEEDLE 3cm 45° ANG, 1cm from TIP	10

SHARP NEEDLE COATED		
P0118	2" (5.08cm) SHARP NEEDLE	12
P0118A	2.5" (6.40cm) SHARP NEEDLE	12

PENCIL		
P0035H	PUSH BUTTON PENCIL w/ COATED BLADE & HOLSTER	25

GROUNDING PAD		
P4500	SINGLE USE REM GROUNDING PAD w/ ATTACHED CORD	12

Item #	Description	Qty.
LAPAROSCOPIC COATED		
P0017	LAPAROSCOPIC STD BLADE 33cm	6
P0017S	LAPAROSCOPIC STD BLADE SPLIT 33cm	6
P0018	LAPAROSCOPIC SPATULA 33cm	6
P0018S	LAPAROSCOPIC SPATULA SPLIT 33cm	6
P0019	LAPAROSCOPIC CURVED BLADE 33cm	6
P0019S	LAPAROSCOPIC CURVED BLADE SPLIT 33cm	6
P0020	LAPAROSCOPIC L-HOOK 33cm	6
P0020S	LAPAROSCOPIC L-HOOK SPLIT 33cm	6
P0020SX	LAPAROSCOPIC L-HOOK SPLIT 44cm	6
P0020X	LAPAROSCOPIC L-HOOK 44cm	6
P0021	LAPAROSCOPIC J-HOOK 33cm	6
P0021S	LAPAROSCOPIC J-HOOK SPLIT 33cm	6
P0021SX	LAPAROSCOPIC J-HOOK 44cm	6
P0021x	LAPAROSCOPIC J-HOOK 45cm	6
P0022	LAPAROSCOPIC NEEDLE 33cm	6
P0022S	LAPAROSCOPIC NEEDLE 33cm	6
P0024	LAPRASCOPIC BALL 33cm	6
P0024S	LAPAROSCOPIC BALL SPLIT 33cm	6
P0026	LAPRASCOPIC O-WIRE 33cm	6
P0027	LAPRASCOPIC THINCURVSPATULA 33cm	6
P0027S	LAPAROSCOPIC THIN CURV SPAT SPLIT 33cm	6
P0100	LAPAROSCOPIC WIRE L-HOOK 33cm	6
P0100S	LAPAROSCOPIC WIRE L-HOOK SPLIT 33cm	6
P0100X	LAPAROSCOPIC WIRE L-HOOK 44cm	6
P0105	LAPAROSCOPIC WIRE J-HOOK 33cm	6
P0105S	LAPAROSCOPIC WIRE J-HOOK SPLIT 33cm	6

LLETZ UNCOATED		
P0410	10mm X 7mm LLETZ SQUARE 13cm	5
P0420	10mm X 5mm LLETZ LOOP 13cm	5
P0430	15mm X 6mm LLETZ LOOP 13cm	5
P0440	20mm X 10mm LLETZ LOOP 13cm	5
P0450	10mm X 10mm LLETZ LOOP 13cm	5
P0460	15mm X 12mm LLETZ LOOP 13cm5	5
P0470	20mm X 12mm LLETZ LOOP 13cm	5
P0480	20mm X 15mm LLETZ LOOP 13cm	5
P0490	20mm X 8mm LLETZ LOOP 13cm	5