

PlumePen® Ultra

Surgical Smoke Evacuation Pencil

FERDINAND^{AP5}
MEDICAL

Telescopic Capture Port with Enclosed Blade

Blade moves in unison with capture port affording versatility in any surgical situation

Integrated design eliminates need for blade replacement saving time and money

Cam Lock System

Prevents capture port movement during surgery

Audible click upon activation/deactivation of locking mechanism

Responsive, Multi-Directional Swivel

Unique, floating ball joint design allows unrestricted hand movement

Integrated Wire/Tubing Management

Ultra-light, PVC free, flexible tubing aids in reducing wrist pull and hand fatigue



"Buffalo Filter has taken surgeon feedback to heart in the new PlumePen® Ultra design bringing surgical smoke evacuation pencils to a new level of function and versatility.

The ability to extend the tip electrode and capture port to almost any desired length while safely securing it in place with the Cam lock system makes for an efficient and versatile device in any surgical situation.

In my opinion, Ultra is the best design to date."

Dr. Dean R. Adams, Colon & Rectal Surgeon
Binghamton, NY



PlumePen® Ultra

Surgical Smoke Evacuation Pencil




7/8 in (22 mm) Cuff
Connects to most surgical smoke evacuators



Cam Lock System
Prevents telescopic capture port from moving


Full Tip Reveal
Unobstructed view of surgical site

Slender, Textured Profile
Improves surgeons grip



Multi-Directional Swivel
Provides for superior hand movement and less wrist drag than competitive products

Telescopic Capture Port
Extends up to 5.3 in (135 mm)



Product Ordering Information

Description	Catalog Number	Quantity
PlumePen® Ultra (Non-Stick ESU Blade) with 3/8 in (9.5 mm) x 10 ft (3 m) hose/tubing & 7/8 in (22mm) filter port adapter	PLPUL2020	20/Case
PlumePen® Ultra (Non-Stick ESU Blade) with 3/8 in (9.5 mm) x 15 ft (4.6 m) hose/tubing & 7/8 in (22 mm) filter port adapter	PLPUL2520	20/Case

Sterile

