

# PlumePen<sup>®</sup> Pro

Surgical Smoke Evacuation Pencil

**FERDINAND<sup>APS</sup>**  
MEDICAL

## Sleek, Slim Profile

*Smaller than competitive surgical smoke evacuation pencils; similar in size to standard ESU pencils*

## Removable Smoke Capture Port

*Allows for easy electrode and capture port replacement*

## Compatible with Standard Electrodes

*Designed for surgeries where longer electrodes are desired*

## Full Blade Tip Reveal

*Provides an unobstructed view of surgical site*

## 360° Swivel

*Allows for fluid hand movement*

## Textured Profile

*Improves surgeon's grip*



Additional capture ports available to accommodate 4 in (102 mm), 5 in (127 mm) and 6 in (152 mm) electrodes.




 **BUFFALO FILTER<sup>®</sup>**  
Filtration Group<sup>®</sup>

# PlumePen® Pro

## Surgical Smoke Evacuation Pencil

The new, improved PlumePen® Pro features an easy to remove x-ray detectable capture port that allows for effortless electrode replacement to accommodate surgeries where longer electrodes are desired. Pro is smaller than competitive surgical smoke evacuation pencils and similar in size to standard electrosurgical pencils.

Pro is a cost effective and suitable alternative for surgeries where a pencil adapter is less optimal.



### Product Ordering Information

Description	Catalog Number	Quantity
PlumePen® Pro (Non-Stick ESU Blade) 10 ft (3 m) tubing	PLPRO4020	20/Case
PlumePen® Pro (Non-Stick ESU Blade) 15 ft (4.6 m) tubing	PLPRO4520	20/Case
Capture Ports (X-ray detectable) (Accommodates 4 in [102 mm] electrodes)	PROTIP420	20/Case
Capture Ports (X-ray detectable) (Accommodates 5 in [127 mm] electrodes)	PROTIP520	20/Case
Capture Ports (X-ray detectable) (Accommodates 6 in [152 mm] electrodes)	PROTIP620	20/Case

### Available Capture Ports



PROTIP420 . PROTIP520 . PROTIP620  
(each sold separately)

