



Provides Precise Spread of Local Anesthetic Solution where it is needed.

Ace Catheter Set

Ace Catheter Set

Ace Catheter Set provides precise spread of local Anesthetic Solution where it is needed.

The Ace Catheter Set is intended for Regional Anesthesia and Peripheral Nerve Block. It is available in various fenestration lengths 2.5cm (8 hole), 6.5cm (15 hole), 10cm (20 hole), 15cm (30 hole) and 30cm (60 hole). Each length of catheters provides precise medication delivery to incisions, subcutaneous, near nerve bundles as well as sub-muscular where it is needed.

- Catheter-19G : 2.5cm (8 hole) / 6.5cm (15 hole) / 10cm (20 hole) / 15cm (30 hole) / 30cm (60 hole)
- Filter - 0.2 μ m
- Introduce Needle - 18G, Tear away sheath
- Ace Lock - Catheter securement
- Syringe - 5cc



- Wide range of fenestration length
- Highly flexible and kink-resistance catheter
- Reinforced with stainless steel coils
- Easy to use introduce needle with tear away sheath
- Gradation on the catheter
- Radiopaque
- Latex free

REF	Specification
TBD-2.5	2.5cm catheter(8 hole), Introduce Needle 18G, 5cc Syringe, Catheter securement, 0.2 μ m Filter
TBD-6.5	6.5cm catheter(15 hole), Introduce Needle 18G, 5cc Syringe, Catheter securement, 0.2 μ m Filter
TBD-10	10cm catheter(20 hole), Introduce Needle 18G, 5cc Syringe, Catheter securement, 0.2 μ m Filter
TBD-15	15cm catheter(30 hole), Introduce Needle 18G, 5cc Syringe, Catheter securement, 0.2 μ m Filter
TBD-30	30cm catheter(60 hole), Introduce Needle 18G, 5cc Syringe, Catheter securement, 0.2 μ m Filter

